

Serstech ChemDash

Installation Guide.

First-time installation guide for ChemDash PC Software Introduction

© Copyright 2024 , Serstech

V1.0

SERSTECH

Contents

1. The	Start	3
1.1.	Setting up accounts for PC using Microsoft Active Directory	5
1.2.	Setting up accounts for PC using Microsoft Azure Active Directory	8
1.3.	Setting up accounts for PC using local accounts	10
1.4.	Instrument users	12
Cre	ating an Arx user with access to Arx Plus	12
1.5.	Setting up accounts for PC using Microsoft Azure Active Directory	14
1.6.	Creating a Serstech instrument user with access to ChemDash 2.0	15

1.7. Creating a Serstech instrument user without ChemDash connect access.



1. THE START

The person Installing the ChemDash application should typically be an IT administrator or equivalent person with administrative rights to the PC. To install the ChemDash software to your personal computer, please follow all steps described below in sequence:

Go to the USB stick and copy-paste the ChemDash file to your local PC.

Setup - ChemDash version 2.0.28	_		\times
Select Destination Location Where should ChemDash be installed?		6	
Setup will install ChemDash into the following folder.			
To continue, click Next. If you would like to select a different folder, click Browse.			
C:\Program Files (x86)\ChemDash	Bro	wse	

At least 684.3 MB of free disk space is required.

Right click on the icon and *Run* the software to go through the installation process.

~ Today (1)			
🖳 ChemDash 2.0.28 setup 210216-1226	2/22/20	Open	102.091 KB
		👽 Run as administrator	
		Troubleshoot compatibility	
		Pin to Start	
		🕂 Scan with Microsoft Defender	
	Ŀ	🖻 Share	
		Open on Mac	
		Show in Finder	
		Share on	>
		Restore previous versions	
		Send to	>
		Cut	
	_	Сору	
		Create shortcut	
		Delete	
		Rename	
		Properties	

178 MB

Once the installation process is finished, double-click on the application's icon to start the application.



4.Then the ChemDash application will open. Now it is time for the IT administrator to set up accounts in the ChemDash application. Depending



on your environment use the most appropriate of the following three possibilities.



1.1. SETTING UP ACCOUNTS FOR PC USING MICROSOFT ACTIVE DIRECTORY

When starting the application for the first time you will automatically belong to the User category SuperAdmin. The SuperAdmin will set up the account for the Administrator of the ChemDash application. This is done by clicking on Admin in the lower part of the menu to the left.

	ChemDash	Welcome noemisocaci						ւհ
≕								
	MEASUREMENTS							
B	REPORTS							
	DEVICES	= 0			Add Hear	Import		
۵	LIBRARIES	· · · ·			Add User	import		ers
		Domain 🔺	📕 Name	. Account	User categ			
۵		Domain	UserName		SuperAd	min		

In the user management access to the ChemDash application is managed through assigning the users to different User categories, with different privileges. The different User categories are explained below:

 SuperAdmin – this category has full rights to all functions of ChemDash The SuperAdmin is typically a person with administrative rights to the computer and IT infrastructure of the company (*i.e.* IT administrator). Thus, this user may not be a user of the ChemDash application apart from installing the first user.

- Admin this category has full rights to all functions of ChemDash special rights: to create a User, to create report templates, to remove and add libraries.
- Active This category has basic access to ChemDash This user has not access to user administration, report templates configuration.
- Disabled This category is an inactive user. An active user that no longer should have access to the system will be put into this category by the Admin. In addition, when creating users on the Serstech instrument without access to ChemDash, these users are assigned user category Disabled.

All users from the active directory are automatically imported, during the start-up, with the User category Disabled. The SuperAdmin's role is to setup a user to administrate ChemDash application. This the done by selecting a user from the list of the imported users by clicking on the username.

4	ChemDash	Welcome Admin					
ĭ.		 ✓ User Management 					
1.	MEASUREMENTS						
	REPORTS						
	DEVICES	Select the device on which to view user prope	erties		-		
	LIBRARIES	System properties					
		Domain: Domain-manual		Name: TestUser			•
		Instrument Properties					
			1				
			USB-drive W				

From the list of available users, the IT administrator will select a person to be Admin for the ChemDash application by clicking on the username. It is possible to have more than one person in the user category Admin.

	ChemDash	Welcome TempAdmin							
≅		< User Management							
	MEASUREMENTS					User Category			
	REPORTS								
	DEVICES	Select the device on which to view user	properties		+				
٩	LIBRARIES					Disabled			
		System properties				Active			
~		Domain: TempAdmin 5-MANUAL		Name: TestUser		Admin			
•						SuperAdmin			
		Instrument Properties							_
					t				
		Modify Scan Settings							
			USB-drive V						
		Change Date And Time	Factory Res						
		Change Language							-
							Cancel	Save	

In addition to setting the User category for the ChemDash user it is also possible to define this user as a user on the Serstech instrument as well. However, this is not necessary, a ChemDash user does not have to have any access to an instrument.

When the IT administrator performing the installation and setting up the first user within the Admin category and is not intended to work with the ChemDash application, it is recommended that this user is deactivated in the application, by the Admin.

The IT administrator's User category is changed from SuperAdmin to disabled. This prevents the user from access to ChemDash. However, the IT administrator can be reactivated as SuperAdmin in the future, when necessary, *e.g.* installation of software updates.

SERSTECH

When starting the application for the first time you will automatically belong to the User category SuperAdmin. The SuperAdmin will set up the account for the Administrator of the ChemDash application. This is done by clicking on Admin in the lower part of the menu to the left.



The SuperAdmin will have to import the users from the active directory into the application by clicking on Import Azure AD users. This will cause the users from the active directory to appear in the list of users in User management.

In the user management access to the ChemDash application is managed by assigning the users to different User categories, with different privileges. The different User categories are explained below;

- SuperAdmin this category has full rights to all functions of ChemDash. The SuperAdmin is typically a person with administrative rights to the computer and IT infrastructure of the company (*i.e.* IT administrator). Thus, this user may not be a user of the ChemDash application apart from installing the first user.
- Admin this category has full rights to all functions of ChemDash special rights: including creating a User, creating report templates, and removing and add libraries.

- Active This category has basic access to ChemDash. This user has no access to user administration, or report templates configuration.
- Disabled This category is an inactive user. An active user who no longer should have access to the system will be put into this category by the Admin. In addition, when creating users on the Serstech instrument without access to ChemDash 2.0 these users are assigned the user category Disabled.

The default User category for the imported users from the active directory is Disabled, thus without access to ChemDash.

From the list of available users, the IT administrator will select a person to be Admin for the application by clicking on the username. From the User category, the Admin is chosen from the drop-down menu. It is possible to have more than one person in the user category Admin. However, it is recommended that further User management is carried out by the Admin of ChemDash.

	ChemDash	Welcome TempAdmin						0	11	սե
ы Ш		< User Management								
	MEASUREMENTS									
B	REPORTS						User Category			
	DEVICES	Select the device on which to view us	er properties	Select Device		Ĵ				
æ	LIRRARIES						Disabled			
	LIDIANILS	System properties					A			
							ACUVE			
		Domain: TempAdmin 5-MANUAL		Name: TestUser			Admin			
							SuperAdmin			
		Instrument Properties								
						ä				
		Privileges								
		Modify Scan Settings	Mount US	B-Drive (read Only)	Account Locked					
		Add Substance	USB-drive	Write Access						
				Application						
		Change Language								
							c	Cancel	Save	

In addition to setting the User category for the ChemDash user, it is also possible to define this user as a user on the Serstech instrument as well. However, this is not necessary, a ChemDash user does not have to have any access to an instrument.

When the IT administrator performing the installation and setting up the first user within the Admin category is not intended to work with the ChemDash application the User Category should be set to Disable, to prevent this user from access t ChemDash.

1.3. SETTING UP ACCOUNTS FOR PC USING LOCAL ACCOUNTS

When starting the application for the first time you will automatically belong to the User category SuperAdmin. All local users available on the PC are automatically imported to the ChemDash application. The local users imported will by default have User category Disabeled and thus no Access to the application. The SuperAdmin will manage the account for the Administrator of the ChemDash application. This is done by clicking on Admin in the lower part of the menu to the left.

	ChemDash	Welcome noemisocaci	Welcome noemisocaci						11	ւհ
≕										
	MEASUREMENTS									
B	REPORTS									
	DEVICES	= 0				Add	ker	Import 4		ers
۵	LIBRARIES									
		Domain 🔺	II Name		# Account		User category			
۵		Domain	UserName				SuperAdmin	1		

In the user management, access to the ChemDash application is managed through assigning the users to different User categories, with different privileges. The different User categories are explained below:

 SuperAdmin – this category has full rights to all functions of ChemDash. The SuperAdmin is typically a person with administrative rights to the computer and IT infrastructure of the company (*f.eg.* the IT administrator). Thus, this user may not be a user of the ChemDash 2.0 application apart from installing the first user.

- Admin this category has full rights to all functions of ChemDash and special rights: to create a User, to create report templates, and to remove and add libraries.
- Active This category has basic access to ChemDash. This user has no access to user administration, or report templates configuration.
- Disabled This category is an inactive user. An active user who no longer should have access to the system will be put into this category by the Admin. In addition, when creating users on the Serstech instrument without access to ChemDash these users are assigned the user category Disabled.

If the SuperAdmin is the only person using the ChemDash application no other users need to be added and the SuperAdmin can remain in this user group. To add additional users the SuperAdmin will manage the users from the list of available by clicking on the username. From the User category drop-down menu the User category is chosen:

	ChemDash	Welcome Admin						
■	MEACHDEMENTC	< User Management						
	REPORTS							
1	DEVICES	Select the device on which to view user prop	erties		Ŧ			
e	LIBRARIES	System properties						
A		Domain: Domain-manual		Name: TestUser		Active		-
Ē		Instrument Properties						
					ā			
		Privileges:						
		Change Date And Time	Factory Rese					

In addition to setting the User category for the ChemDash 2.0 user, it is also possible to define this user as a user on the Serstech instrument as well. However, this is not necessary, a ChemDash user does not have to have any access to an instrument.

1.4. INSTRUMENT USERS

Creating an Arx user with access to Arx Plus.

To be able to create users on a Serstech instrument, the device must have been connected to ChemDash application.

	ChemDash	Welcome noemisocaci			
×	MEACHDEMENTS	ζ User Management			
-	MEASUREMENTS				
E	REPORTS				
	DEVICES	Select the device on which to view user properties Select Device			
۵	LIBRARIES	System properties			
		Domain: N Domain -MANUAL Name: User User Admin			-
6					
		Instrument Properties			
		PIN-code PIN Code Expiry:			
		Privileges: Account:			
		Modify Scan Settings Mount USB-Drive (read Only)			
		Add Substance USB-drive Write Access			
		Change Date And Time Factory Reset			
		Change Language			
		Canc	el	Save	

Click on the user from the list of users in the User management to be able to edit.

Select the user category, Active or admin, from the drop-down list.

Select the instrument from the drop-down list in Select Device.

Configure instrument.

Set PIN-code, to be used for Serstech instrument log-in.





Set PIN-code expiry date and time.



Set privileges for the user on the instrument by ticking the options in the list of privileges.

Privileges:		Account:
Modify Scan Settings		
Add Substance	USB-drive Write Access	
Remove Measurements	Pair With Application	
Change Date And Time		
Change Language		

Save.

The user is sent to the instrument in the Devices section.

Note that the ChemDash Admin can remove access to the application by changing the User category to Disabled and thus preventing log on to ChemDash.

1.5. SETTING UP ACCOUNTS FOR PC USING MICROSOFT AZURE ACTIVE DIRECTORY

If you have an account using Active Directory, you need your company's IT administrator to import you into the application the first time. The IT administrator is an admin outside of the system who when installing the ChemDash software will have the access privilege SuperAdmin. As a SuperAdmin the IT administrator will import all available users in the Azure active directory by clicking on the button Import Azure AD users.

	ChemDash	Welcome noemisocaci				©	1	ılı
≕								
03	MEASUREMENTS							
B	REPORTS							
	DEVICES							
۵	LIBRARIES	÷ < •			Add User			ers
		Domain 🔺	🔛 Name	Account	User category			
ê		Domain	UserName		SuperAdmi	n		

From the list of available users, the IT administrator will select a person to be Admin for the Arx application by clicking on the username. It is possible to have more than one person in the user category Admin.

4	ChemDash	Welcome Admin							
×	MEACUDEMENTS	< User Management							
<u> </u>	MEASUREMENTS								
	REPORTS								
	DEVICES	Select the device on which to view user properties				-			
e	LIBRARIES	System properties							
		Domain: Domain-manual		Name: TestUser					-
ê		Instrument Properties							
						Ċ,			
		Privileges:							
		- Modify Scan Settings		Drive (read Only)					

In addition to setting the User category for the ChemDash user it is also possible to define this user as a user on the Serstech instrument as well.

SERSTECH

However, this is not necessary, a ChemDash user does not have to have any access to an instrument.

When the IT administrator performing the installation and setting up the first user within the Admin category is not intended to work with the ChemDash application, the User Category should be set to Disable, to prevent this user from access to ChemDash.

1.6. CREATING A SERSTECH INSTRUMENT USER WITH ACCESS TO CHEMDASH 2.0.

To be able to create users on a Serstech instrument the instrument must have been connected to ChemDash application.

	ChemDash	Welcome noemisocaci		
×	MEACHDEMENTC	∠ User Management		
-	MEASUREMENTS			
1	REPORTS			
	DEVICES	Select the device on which to view user properties Select Device -		
۶	LIBRARIES			
		System properties		
		Domain: N Domain -MANUAL Name: User User Admin		-
ê				
		Instrument Properties		
		PIN-code:		
		Privileges: Account:		
		Mount USB-Drive (read Only)		
		Add Substance USB-drive Write Access		
		Remove Measurements Pair With Application		
		Change Date And Time Factory Reset		
		Cance	Save	

Click on the user from the list of users in the User management to be able to edit.

Select user category, Active or admin, from the drop-down list.

Select instrument from the drop-down list in Select device.

Configure instrument.



Set PIN-code, to be used for Serstech instrument log on.

Add PIN-code					×
	3	4	4	0	
	Cancel		Ad	id PIN-code	

Set PIN-code expiry date and time.



Set privileges for user on the instrument by ticking the options in the list of privileges.

Privileges:		Account:
Modify Scan Settings		
Add Substance	USB-drive Write Access	
Remove Measurements	Pair With Application	
Change Date And Time		
Change Language		

Save

The user is sent to the instrument in the Devices section.

Note that the ChemDash Admin can remove access to the application by changing the User category to Disabled and thus preventing logging on to ChemDash.

1.7. CREATING A SERSTECH INSTRUMENT USER WITHOUT CHEMDASH CONNECT ACCESS.

In the User management, a new user is added through Add User.

Type a name and select a User category. Disabled is an appropriate user category since this user will not have access to ChemDash.

Add User			×
Name: New User	Disabled	-	
		Add User	1

Once added the user will appear in the list of users in the User management. The family name in the list will be the name of the local domain hyphened with the manual.

Configure the created user for an instrument and select the user by clicking on the name in the list displayed in the User management.



ChemDash	Welcome noemisocaci			
MEACUDEMENTE	ζ User Management			
MEASUREMENTS				
REPORTS				
DEVICES	Select the device on which to view user properties Select Device +			
LIBRARIES				
	System properties			
	Domain -MANUAI Name User User			÷
	Instrument Properties			
	PIN-code:			
	Privileges Account:			
	Moufify Scan Settings Mount US8-Drive (read Only) Account Locked			
	Add Substance USB-drive Write Access			
	Change Language	_		1
	Cance	ł	Save	
	ChemDash MEASUREMENTS REPORTS DEVICES LIBRARIES	ChemDash Velcome nomisocai MEASUREMENTS RPORTS DEVICTS LIBRARIS DEVICTS Domain Domain MANUAL Name User User Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument Properties Instrument	ChemDash Wetome neemisocasi ALMAN MEASURIGHARINS REFORTS DEVICES DEVICES DEVICES Demain Domain MANULL Name User User Admin Instrument Properties Pint-code Pint-code <th>ChemDash Wetcome consistanci ALEXAUREMENTS REFORTS DEVICES LIBRARIES System properties System properties System properties Imatriament Properties Miniter Standard Standard Pindeges Individual Standard Ind</th>	ChemDash Wetcome consistanci ALEXAUREMENTS REFORTS DEVICES LIBRARIES System properties System properties System properties Imatriament Properties Miniter Standard Standard Pindeges Individual Standard Ind

Select the user category Disabled from the drop-down list.

Select the instrument from the drop-down list in Select Device.

Configure instrument.

Set PIN-code, to be used for Serstech instrument log-on.



Set PIN-code expiry date and time.





Set privileges for the user on the instrument by ticking the options i the list of privileges.

Privileges:		Account:
Modify Scan Settings	Mount USB-Drive (read Only)	Account Locked
Add Substance	USB-drive Write Access	
Remove Measurements	Pair With Application	
Change Date And Time	Factory Reset	
Change Language		

Save

The user is sent to the instrument in the Device section.