

SERSTECH CHEMDASH 2.0 Manage - Analyze - Share


Serstech Chemdash 2.0 is a powerful software tool that manages and connects the Serstech Raman devices with the rest of your organization.

Identification and quick decision making are often the main focus when working with possibly hazardous or illegal substances. Having fast and reliable instruments in place is important for the field operator. Serstech ChemDash 2.0 is the software that makes the identification solution complete for all roles in the team.

Manage devices, users, measurements and reports from one single software. The devices can be prepared with users and other settings, data can be collected and analyzed, and the results can be shared with the rest of the organization. The ChemDash administrator can create and manage their own libraries, making it possible to tailor sets of substances, based on the measurements. Complete the solution for substance identification and decision making, all through ChemDash.

ChemDash is a powerful tool for all the above and still it only takes seconds to learn. ChemDash 2.0 Lite is a free version made for the small organization, where one instrument can be connected to the system. ChemDash 2.0 Pro brings more possibilities in terms of system and device customizations and allows several devices to be connected and managed.

"ChemDash is a powerful tool and still it only takes seconds to learn."



New user interface makes navigation intuitive and efficient. Create libraries, manage measurements, update user settings and share your results with the rest of your organization.



For questions regarding the Serstech Partner Program, visit our website www.serstech.com/contact or contact one of our regional sales managers.

SERSTECH

ChemDash	God morning, Stefan!	0	□ ⑦ 11.	
≓< ★ measurements	< Polyetylen 19 (<u>edit)</u>	Spectrum		
REPORTS REPORTS REPORTS DEVICES USERS USERS ADMIN	Polycybien 19 (ggl) Substance basic information Substance basic information			• 17550
Measuremen	t information including spectrum w	ith zoom functionality.	00 800 10	00 1200 1400 1600

ChemDash	God morning, Stefan!	ار © 💵 🧿 ۱	ı.				
₩ MEASUREMENTS	X Title Time Device Scan1D Type Q III						
	📄 🔹 📩 🗄 Title 🔺 🗄 Spectrum 🔛 Time 🔛 Device						
UBRARIES	Alprozolam, 194 – 🥠 2020 02 24, 11:09 Device 1	355 Narcotic **					
LUSERS	Amorbabital, 195 2020 02 24, 11.09 Device 1	356 Not regulated ••					
admin	Amphetamine sulfate, 196 JML 2020 02 24, 11:09 Device 1	357 NPS **					
	Barbital, 197						
	Bromazepane, 198						
	Bromazepane, 199						
	Bromazepane, 200 - M. 2020 02 24, 11:09 Device 2						
	A. 5 5 6 1 4 100-199 of 1547 > 3 4	4 5 >>> Show	•				

Easy to manage measurements and devices. Analyze the data and share with the rest of your organization.

Information	ChemDash 2.0 Lite	ChemDash 2.0 Pro
Connectivity	USB	USB, Wi-Fi
Create reports	Yes	Yes
Customize reports	-	Yes
User management	Basic	Advanced
Number of users	1	Up to 100
Device management	Yes	Yes
Customize device settings	-	Yes
Maximum number of devices	1	Multiple