

Chemdash 2.7.23

This release includes:

- Inclusion of noisy measurements in libraries is no longer permissible
- Exported libraries no longer feature password protection; consequently, the json file is accessible without authentication
- The requirement to periodically reconnect to Active Directory has been eliminated
- Downloading of old measurements is now possible
- The device network tab has been removed
- Exported logs are ordered chronologically
- Downloaded measurements are presently arranged based on scan time, with the most recent being prioritised for download
- Better and more responsive spectra graph
- Selection of measurement images is restricted exclusively to image files
- Improved system user detection method
- Download measurements progress bar is significantly improved
- Maximum of 25 substances that can be added into a report
- Maximum of 10 substances that can be compared at once
- Exporting is limited to published libraries
- Animations for report loading and downloads have been incorporated
- Censored Wi-Fi password for increased security
- Buttons with lacking or exceeded clickable surface has been resolved
- Libraries lacking an English name can now be exported
- Modifying any setting will no longer result in a deduction of PIN code expiry time
- The language selector is now displayed within the device system settings
- During measurements download, the PIN code field will no longer reset to the current value after a certain period

This constitutes a summary of various enhancements and bug fixes for Chemdash. If you have specific inquiries, please don't hesitate to reach out to our support team.