

# Firmware 6.8.0

## This release includes:

- Enhanced fluorescence removal
- Corrected certain peak-detection cases where an obviously good spectra could still result in a 'No match'
- Created labels are now promptly displayed on the result and scan history screens
- The expiration time for PIN code is now consistently enforced
- User privileges are now configured accurately
- Factory reset function is operational when the instrument is connected to Chemdash
- Added substances to the Narcotics library:  
Promethazine, BOC-ANAM, Primidone, Acamprosate calcium, 1-benzyl-piperazine
- French translation has been integrated for the entire Demo library

This constitutes a summary of various enhancements and bug fixes for Arx and i100.  
If you have specific inquiries, please don't hesitate to reach out to our support team.